
DENTAL SINTERING EQUIPMENT

F10+

Zirconia Sintering Furnace

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QUICK GUIDE

Each language section includes the same technical content. Use the table of contents at the beginning of your chosen section. Sintering programs · installation · operation · maintenance.

F10+

Zirconia Sintering Furnace — Instruction Manual

ENGLISH

EN — TABLE OF CONTENTS

Overview of the sections and principal procedures in this manual.

1. Overview	2. Technical Parameters
3. Parts List	4. Introduction to Equipment Structure
5. Equipment Installation	6. Startup
7. Equipment Operation 7.1 General Steps (Running a Program) 7.2 Program Curve Setting Steps	8. Product Highlights
9. Notes	10. Maintenance

1. Overview

Thank you very much for purchasing and using Yucera's products. This series of sintering furnaces adopts an ergonomic appearance design, with a novel design and reasonable structure. Utilizing modern mechanical manufacturing and processing equipment, the furnace is produced with exquisite workmanship and superior quality.

The high-temperature furnace liner is made of imported alumina lightweight fiber material with excellent insulation performance, environmental protection, and pollution-free performance. The control interface is a 7-inch high-definition true-color LCD touch screen, with graphic display and easy operation. The imported PID intelligent temperature control offers high temperature-control accuracy and reliable operation.

Please read this operation manual carefully before using this product.

2. Technical Parameters

Specification	Parameters
Temperature control accuracy	$\pm 1^{\circ}\text{C}$
Touch-screen size	7" true color
Rated power	$\leq 2 \text{ kW}$
Number of programs	24 programs
Heating rods	4 pcs, 1850-type high-purity silicon-molybdenum rods
Sensor type	S-type, double platinum-rhodium
Maximum temperature	1560°C
Operating temperature	$\leq 1550^{\circ}\text{C}$
Heating rate	0–30 $^{\circ}\text{C}/\text{min}$ (0–500 $^{\circ}\text{C}$) 0–10 $^{\circ}\text{C}/\text{min}$ (500–1200 $^{\circ}\text{C}$) 0–12 $^{\circ}\text{C}/\text{min}$ (1200–1530 $^{\circ}\text{C}$)
Furnace size	$\varnothing 110 \text{ mm} \times 130 \text{ mm}$
Overall dimensions	Width 37 × Depth 53 × Height 78 cm
Equipment weight	Approx. 45 kg

3. Parts List

No.	Name	Unit	Qty
1	Zirconia sintering furnace	set	1
2	Zirconia slow-burning crucible	pcs	2
3	Special fast-burning crucible	pcs	1
4	Crucible tongs / crucible tongs for fast-burning crucible	pcs	each 1 pc
5	Sintered beads	bag	1
6	Stove (base plate)	pcs	1

4. Introduction to Equipment Structure

F10+ furnace structure

① Sintering furnace body ② Bottom tray ③ Touch-screen control panel

5. Equipment Installation

Installation tips:

- Open the packing wooden box, check whether the equipment is in good condition, and check whether the accompanying accessories are complete according to the packing list.
- Place the furnace body on a dry, ventilated, and dust-free indoor platform. The distance from surrounding objects should not be less than 25 cm, so the equipment has good ventilation and heat-dissipation conditions during normal operation.

- Pay attention to the heat resistance of equipment placed nearby. The device releases a certain amount of heat to the surrounding environment during use; although the heat value is within a safe range, it is possible that peripheral equipment casings could discolor.
- It is strictly forbidden to place flammable or explosive materials around the furnace body.

When installing the equipment's power cord, pay attention to the power-supply mark on the cord and do not use a different supply voltage than indicated (domestic supply is 220 V / 50/60 Hz), to avoid damage to the control system. Turn off the power supply when not in use.

Note: This device is a 2 kW high-power device. Please plug it into a socket of guaranteed quality, and ensure the socket's power line has a cross-section greater than 2.5 mm².

6. Startup

The main power switch is located on the side of the machine screen — push it up to turn the equipment on.

7. Equipment Operation

A. The touch-screen startup interface is as follows:

Startup screen — tap "Log in" to open the main display

Main display — program, temperature curve, running state, PV/SV, and control buttons (Rise, Pause, Decline, Start, Stop, Menu)

7.1 General Steps

1. Select the curve and confirm the curve

Program selection window (P1–P24)

- P1–P16 on this machine are slow-burning curves; the starting temperature for platform lowering is $\leq 500^{\circ}\text{C}$.
- P17–P24 are special curves for fast burning; the starting temperature of the stove base is $\leq 1000^{\circ}\text{C}$. Be sure to use the special crucible for fast burning.
- P17–P20 are fast-burning reference curves pre-set for customers to choose from.

2. Check that the crucible is centered on the stove top

3. Click "Start", then click "OK" to begin — this enables one-key start and run

"Please confirm the run" dialog — tap "Determine" to start


7.2 Program Curve Setting Steps

1. Tap "Menu", then tap "Curve setting"


Main menu — Main control interface / Curve setting / Set / Operation record / Historical curve / Data export / Alarm record / Exit / User switching

2. Program-settings page


Curve-settings table (end temperature, time, rate per segment) and process selection panel

3. Customers can rename the curve themselves

"Please enter name" dialog

On-screen keyboard for entering the new curve name

8. Product Highlights

Yucera's zirconia fast-slow integrated sintering furnace is a new-generation presentation. The heat-preservation furnace, made of imported raw materials, is equipped with 1850-type high-standard silicon-molybdenum rods, giving every crown excellent performance while using only a maximum power of 2 kilowatts — about 2.3 kWh. Sintering of about 15 crowns can be achieved in as fast as 90 minutes.

The traditional 10-hour slow-burning cycle can be performed by stacking two layers of 100 mm-diameter crucibles, and the sintering can be completed using only 8.5 kWh. It can run fast or slow, in the factory or chairside, saving energy and electricity — and the whole machine is guaranteed for one year.

1. Integrated fast-burning and slow-burning furnace, which reduces the equipment-purchase cost for the laboratory.
2. Sintering of a front-tooth single crown can be completed in as little as 90 minutes, ready for polishing in 2 hours. The full sintering effect can be achieved within 5–7 hours in daily use, with a transparent, deformation-resistant color. Even sintering carried out late at night will not affect the next day's grinding schedule, improving the factory's production efficiency.
3. This furnace offers 24 sintering programs, 4 of which are preset curve programs provided by the furnace manufacturer for customer reference.
4. A laboratory processing a large volume of restorations can sinter three times a day without affecting normal shipments; color remains accurate and pollution-free.

9. Notes

- [A] When the furnace is used for the first time, or after a long period of non-use, it should be fired following the normal procedure. The furnace temperature must not exceed the rated temperature, so as not to damage the heating element and furnace lining.
- [B] After the program finishes running, it will automatically enter the stopped state.
- [C] When the equipment is not in use, the lifting platform should be raised to the top to keep the furnace cavity dry.

10. Maintenance

- [A] During high-temperature sintering, do not open the furnace. This is not only dangerous, but the heating element is also easily bent and deformed by sudden temperature changes, and the furnace itself is prone to cracking and damage.
- [B] After sintering, confirm whether the temperature inside the furnace is below 300°C. If it is above 300°C, do not cut off the power supply — this would stop the cooling fan, and the furnace's surface temperature could rise sharply, damaging the furnace body.
- [C] When the furnace is used for the first time, or after a long period of non-use, it should be baked at about 120°C for 1 hour, then at about 300°C for 2 hours, before use, to prevent the furnace from cracking. The furnace temperature must never exceed the rated temperature, so as not to damage the heating element and furnace lining. It is forbidden to pour liquids or molten metals directly into the furnace — keep the furnace clean.
- [D] The furnace body uses silicon-molybdenum rods as heating elements. Due to their physical characteristics, these rods are very brittle at room temperature. Therefore, the furnace body must not be disassembled or moved carelessly once the heating elements are installed.
- [E] When using a cold furnace, since the furnace is cold it must absorb a large amount of heat, so the temperature-rise rate of the low-temperature section should not be too fast, and the rate difference between temperature sections should not be too large. When setting the heating rate, fully consider the physical and chemical properties of the sintered material, to avoid spattering and contamination of the furnace.
- [F] Regularly check whether the contacts of the temperature-control system's electrical connections are in good condition, paying special attention to the connection and tightness of each connection point of the heating element.
- [G] When silicon-molybdenum rods are used as heating elements, it is not suitable to operate for long periods within the 400–700°C range, as the rods will oxidize at low temperature under those conditions.

The furnace is suitable for the following working conditions:

- Ambient temperature between -10°C and 75°C.
- Relative humidity of the surrounding environment not exceeding 85%.

- No conductive dust, explosive gas, or corrosive gas that could seriously damage metal or insulating materials around the furnace body.
- No noticeable tilt, vibration, or bumping.

Warranty

Provided the user complies with the storage, use, installation and transportation regulations, if the product is damaged due to a product-quality issue and cannot work normally within 12 months of the delivery date, our company will provide the user with warranty service for the whole machine (man-made damage excluded). After the warranty period expires, our company will continue to offer paid lifetime maintenance according to user requirements.

Note: the heating element is a consumable item, but the company still guarantees it for 12 months. If you encounter any problems during use, please contact us promptly.

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