## DENTAL MILLING EQUIPMENT

# YRC-5X

## Five-Axis Dental Milling Machine

Operation Manual · Manual de Operação · Manual de Operación

EN · PT-BR · ES

## MULTILINGUAL EDITION

YRC-5X · Five-Axis Dental Milling Machine

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### QUICK GUIDE

Each language section includes the same technical content. Use the table of contents at the beginning of your chosen section. Installation · machining operation · tools · maintenance.

# YRC-5X

Five-Axis Dental Milling Machine — Operation Manual

---

ENGLISH

---

# EN — TABLE OF CONTENTS

Overview of the sections and principal procedures in this manual.

## 1. YRC-5X Milling Machine

- 1.1 Equipment Introduction
- 1.2 Host Operation Panel
- 1.3 Working Cabin Door
- 1.4 Working Cabin
- 1.5 Suction and Blowing System
- 1.6 Technical Parameters

## 2. Equipment Transportation and Installation

- 2.1 Handling, Packaging and Storage
- 2.2 Parts List
- 2.3 Installation Requirements
- 2.4 Device Installation
- 2.5 Removing the Working-Chamber Filling
- 2.6 Pneumatic Unit Installation
- 2.7 Installing the Vacuum System
- 2.8 Connecting the Circuit

## 3. Operating the Equipment

- 3.1 Preparation Before Operation
- 3.2 Powering On the Equipment
- 3.3 Equipment Testing
- 3.4 Overview of the Milling Process
- 3.5 Inserting / Replacing the Burs
- 3.6 Burs Information
- 3.7 Installing and Removing the Fixture
- 3.8 Program Interruption and Stop

## 4. Operating System Interface


- 4.1 Status Display
- 4.2 Button Functions
- 4.3 Coordinate Display Area
- 4.4 Tool Position

## 5. Equipment Maintenance and Cleaning

- 5.1 Safety Notes
- 5.2 Accessory Replacement
- 5.3 Routine Cleaning
- 5.4 Material Fixture Cleaning
- 5.5 Cleaning the Tool Library and Tools
- 5.6 Cleaning the Spindle Clamp
- 5.7 Checking Hoses, Cables and Connections
- 5.8 Checking the Pressure-Regulating Filter
- 5.9 Maintenance of Screw and Guide Rail
- 5.10 Maintenance Schedule
- 5.11 Equipment Calibration

# 1. YRC-5X Milling Machine


## 1.1 Equipment Introduction

*Schematic diagram of the YRC-5X*

1. Cutting machine 2. Cooling fan 3. Cabin door 4. Observation window of the working cabin 5. Digital display barometer / USB interface 6. Host system 7. Power on 8. Emergency stop 9. Processing alarm display bar


The YRC-5X five-axis milling machine is designed for the convenient realization of digital processing in the field of dental technology, and can process different types of materials.

Connection board on the back of the device

1. Power inlet connector 2. Network port 3. Dust suction port 4. Air intake hole 5. Vacuum cleaner external interface 6. Double-cup air pressure filter

## 1.2 Host Operation Panel

Operation panel

1. System display 2. Power on/off 3. System emergency stop

## 1.3 Working Cabin Door

The working-cabin door must be closed while the machine is running, to prevent personal injury.

*Note: Beware of pinching when opening or closing the hatch.*


- When opening the hatch, hold the hatch handle with your hand and pull it outward.
- To close the hatch, place your hand on the lower edge of the hatch and push upward.

YRC-5X opening-door rendering

## 1.4 Working Cabin

Compartment for processing materials.

YRC-5X work cabin

1. Tool-setting device 2. Tool magazine 3. Blank fixture 4. Four rotating shafts 5. Five rotating shafts 6. Spindle

## 1.5 Suction and Blowing System

The suction and blowing system includes:

- Gas blown from the spindle and blow-pipe.
- External extraction system (suction device, suction hose).

During normal operation, gas is continuously blown into the working chamber, and the connected dust-suction device draws the debris and dust out of the working chamber through the dust-suction hole. A good extraction and blowing system prevents damage and wear to sensitive parts of the machine.

*Note: The suction and blowing system cannot replace regular machine cleaning. If the machine is not cleaned regularly, its service life will be significantly reduced.*

## 1.6 Technical Parameters

| Specification          | Parameters                                           |
|------------------------|------------------------------------------------------|
| Number of axes         | Five-axis linkage                                    |
| Cutting method         | Dry milling                                          |
| Cutting material       | Zirconia block, PMMA, wax, resin, composite material |
| Block size             | Straight warp 98 mm, thickness 10–30 mm, disc        |
| Stroke range           | X/Y/Z: 200-120-120 A: 360° B: +35°, -92°             |
| Spindle                | 0–60000 rpm / 500 W                                  |
| Drive                  | Full servo motor                                     |
| Tool (bur)             | Quantity: 7 Diameter: 4 mm                           |
| Dimensions             | Left/right: 530 Front/back: 650 Up/down: 760 (mm)    |
| Weight                 | Total mass of the cutting machine: 145 kg            |
| Required air pressure  | 6.0–6.5 bar                                          |
| Compressed air flow    | Approx. 50 L/min                                     |
| Minimum suction volume | Approx. 3000 L/min                                   |
| Power                  | 800 W                                                |
| Power supply           | 220–230 V AC, 50/60 Hz                               |

## 2. Equipment Transportation and Installation

---

### 2.1 Handling, Packaging and Storage

*Warning: In order to ensure the normal use of the equipment in the future, operators and maintenance personnel are required to be competent in handling the transport packaging — please read the following precautions carefully.*

#### 2.1.1 Handling Safety Rules

- Pay attention to the center of gravity of the equipment. Do not shift the center of gravity during transportation, otherwise the device may tilt or fall.
- Pay attention to the weight and direction marked on the outer packing box, and try to keep balance during transportation.
- When using hoisting machinery, always keep the mechanical force around the center of gravity, and minimize shaking.
- When lifting and lowering, always watch for tilting, and reduce packaging height as much as possible if conditions permit.
- Only use the designated stress area during the moving process.
- The outer packaging material may only be removed just before the equipment is installed.

#### 2.1.2 Packing

The transportation package is designed according to the expected transportation method and conditions, and is specially designed to prevent physical damage during transportation to the greatest extent possible. Therefore, it may only be removed just before installation. When removing the transport package, try to preserve its integrity in case it needs to be returned to the factory for reuse.

#### 2.1.3 Logo on the Outer Packaging

- The equipment is a high-precision instrument with many vulnerable internal parts — handle it with care during transportation.
- There are many electronic components and mechanical parts inside that need to be kept away from moisture.
- When temporarily placing or storing the equipment, place it in the direction of the arrow.
- Do not let the box tilt or roll over during transportation and handling.

## 2.2 Parts List

First check the delivery list to see if it is complete. Unpack the machine and check the following items:

| No. | Product name                 | Specification                                                                | Unit | Qty |
|-----|------------------------------|------------------------------------------------------------------------------|------|-----|
| 1   | YRC-5X                       | 530-650-760                                                                  | set  | 1   |
| 2   | Dust suction air pipe        | 2 m                                                                          | pc   | 1   |
| 3   | 8 mm air pipe                | 5 m                                                                          | pc   | 1   |
| 4   | TY tee                       | 8-turn double 8 /<br>6-turn double 8                                         | pc   | 1   |
| 5   | Hoop                         | Vacuum pipe hoop                                                             | pc   | 1   |
| 6   | Power cable                  | 3.5 m                                                                        | pc   | 1   |
| 7   | Screwdriver                  | Phillips and flat-blade                                                      | set  | 1   |
| 8   | Allen wrench                 | 1–10                                                                         | set  | 1   |
| 9   | Tool burs                    | ZIR 2.0 / ZIR 1.0 / ZIR<br>0.6 / 1.5R / PMMA 2.0<br>/ PMMA 1.0 / PMMA<br>0.6 | pc   | 7   |
| 10  | Pressure regulator<br>filter | AFR2000 (mounted on<br>the machine)                                          | pc   | 1   |
| 11  | Wax disc                     | 98 mm                                                                        | pc   | 1   |
| 12  | Lifting handle               | 20–160                                                                       | pc   | 1   |


## 2.3 Installation Requirements

The installation site must meet the following conditions:

- The installation site must be firm, with a flat surface able to bear the weight of the equipment (approx. 150 kg).
- Ideal indoor temperature: 18°C–25°C, maximum 32°C.
- The installation location must be dust-free.
- Relative air humidity up to 80%.
- AC power supply of 220–240 V, 50/60 Hz.
- Compressed-air supply according to machine requirements: 4.5–7.5 bar.
- Reserve enough space around the installation: 10 cm on top, left, right and back of the machine.

*Note: Plug the unit's power cord into a separate fused circuit, or make sure no other equipment shares the same power source, to prevent voltage fluctuations. Strong voltage fluctuations can interfere with control and cause system failure.*

## 2.4 Device Installation

Device connection diagram — Power supply / Cable / Vacuum / Air intake / Electric vacuum

## 2.5 Removing the Working-Chamber Filling

Before operating the device for the first time, the filling in the working chamber must be removed. This filling protects the spindle and machine structure during transport.

- Remove the shipping stuffing from the working compartment.
- Clean up any remaining stuffing debris in the working compartment.
- Keep the filling safe for later transportation and use.

## 2.6 Pneumatic Unit Installation


- During installation and servicing, ensure no compressed air passes through the regulator filter.
- After installing the pneumatic hose, before venting, check the hose is securely inserted into the correct connector and undamaged.

*Note: The compressed-air input must be oil-free and dry. Contaminated compressed air (containing water, oil, etc.) may damage the spindle and electrical components.*

Connect the machine to a compressed-air supply only through the supplied pressure-regulating filter.

### 2.6.1 Introduction to the Pressure-Regulating Filter

The compressed air is connected to the machine through a pressure-regulating filter, which adjusts the pressure input to the machine and partially separates moisture from the gas.

*Regulator filter*

1. Pressure adjustment knob 2. Fixed seat 3. External compressed-air interface 4. Pressure gauge 5. Water separator 6. Drainer 7. Machine connection interface

## 2.7 Installing the Vacuum System

### 2.7.1 Vacuum Cleaner Requirements

Only use vacuum cleaners with the following characteristics:

- Specially designed for dentistry.
- Equipped with a class-M filter.
- Suited to the working environment of the machine tool.
- Fitted with a mounting device to prevent operator injury.
- Minimum displacement of 3000 L/min.

### 2.7.2 Vacuum Cleaner Installation

Before installing the vacuum cleaner, read its installation and operation instructions carefully.

- Check whether the vacuum cleaner's suction port matches the bellows (inner diameter 49 mm); if not, use an adapter.
- Insert the bellows into the suction port of the vacuum cleaner and make sure the connection is firm.
- Insert the bellows into the suction port of the machine and make sure the connection is firm.
- If you want the machine to automatically switch the vacuum cleaner on and off, connect it to the machine through a linkage line.

## 2.8 Connecting the Circuit

- Connect and lock the aviation connector on the cutting machine and the corresponding aviation connector on the main unit, in sequence, using the connecting cable.
- Plug the supplied power cord into the power connector on the console's connection panel.
- Plug the power cord into a grounded, overload-protected outlet.

*Caution: Severe voltage fluctuations can damage the device and interrupt the control unit, possibly causing system failure. Plug the device's power cord into a dedicated circuit, or ensure nothing else that could cause large voltage fluctuations is connected when powered on. If severe fluctuations cannot be avoided, install a voltage stabilizer.*

## 3. Operating the Equipment

### 3.1 Preparation Before Operation

- Ensure the equipment placement surface is horizontal and stable, not prone to vibration or shaking.
- Ensure all components are properly installed and connected.
- Ensure all spare parts are firmly fixed.
- Ensure all control systems are installed and operational.

### 3.2 Powering On the Equipment

- Ensure the equipment is fully installed, with sufficient compressed air supplied, and close the working door.
- Turn on the power switch on the host backplane.
- Press the switch button on the host operation panel to start the equipment.
- Home the equipment — after the control system starts fully, the equipment will begin homing automatically.

*Note: The working-cabin door must be closed whenever the equipment is started and running. Only after the homing action is complete and the reminder window on the display is closed can the next step be performed.*

### 3.3 Equipment Testing

After first installation or reinstallation (after transportation), the basic performance of the equipment must be tested.

- Put the burs into the machine.
- Without installing a zirconia block or wax disc on the fixture, load any program and run the equipment.
- Observe whether the equipment runs normally (tool changing, blowing, movement of each axis, etc.).

- Confirm the equipment has no problems and can process normally. If it cannot operate normally, contact the manufacturer.

### 3.4 Overview of the Milling Process

- Before starting processing, read the operation manual and safety instructions carefully.
- Check the burs and fixture thoroughly for damage before operating. Damaged burs or fixtures may damage the equipment.
- The processing program is an "nc" file, typeset via CAM software such as hyperDENT. This file contains all the information required for processing. Copy the program into the operating system via USB or via the network.


Follow this workflow with the YRC-5X:

1. Select the corresponding block size and material type in the hyperDENT software for typesetting.
2. Insert the burs into the equipment.
3. Put the material block into the equipment and clamp it tightly.
4. Click the "Load" button to select the program to run from the program list.
5. Click "Start-up" to begin processing.

### 3.5 Inserting / Replacing the Burs

- Incorrect burs will damage the spindle's collet or tool magazine.
- The bur needs a sufficiently large chamfer on the shank, otherwise the collet is easily damaged.
- Do not use burs with a blade diameter larger than the shank diameter ( $\varnothing 4$  mm).

You can insert 7 burs into the tool changer as shown. The equipment can automatically change the burs during processing.


Tool magazine — bur positions 1–7

You can install the burs in two ways:

**Manually insert the burs into the tool magazine:**

1. Make the spindle stop at a safe position, enter the second page of the operating system, pinch the middle part of the bur on the spindle with your hands, and click "Loose/Clamp" to take it out.
2. Insert the new bur into the spindle, making sure the bur-fixing copper ring is tightly attached to the spindle chuck, then click "Loose/Clamp" to clamp the bur.
3. Click "Auto T" to detect the length of the current new tool and automatically record it into the system.

*Bur shank — milling / copper ring / clamping*

*Note: After the clamping tool is released every time — whether or not the bur is replaced — the tool length must be checked, otherwise the workpiece may be damaged or the tool may break.*


#### Through the tool magazine:

1. Take out the old burs one by one from the tool holder.
2. Put the new bur into the corresponding tool seat correctly according to the serial number; the positioning copper ring should be completely embedded in the tool seat.
3. Confirm the correct placement order and position of the burs.

*Note: Do not hit the bur during placement to prevent damage.*

- When the running time exceeds the preset time, or the workpiece is damaged, replace the bur promptly.
- The bur should be replaced in the tool magazine whenever possible, unless the bur on the spindle is damaged and its length cannot be detected.
- Bur replacement must be performed by trained personnel familiar with the equipment and software.
- The bur must be placed in the corresponding tool holder according to the preset specifications.
- The current default tool number is shown on the operating system. If all burs are returned normally to the tool magazine, the current tool number is displayed as "0".

A length check is performed every time a bur is removed or before it is put back in. The length deviation cannot exceed the preset range; if no bur is detected or the deviation exceeds the range, the system displays an alarm.

*System alarm message screen*

If a tool error occurs during processing, the program must be restarted after checking or replacing the bur.

### 3.6 Burs Information

The table below shows the bur specifications provided by the manufacturer for the YRC-5X.

*Note: Using non-manufacturer burs may lead to abnormal operation or damage to the working parts.*

| Tool No. | Specification                   |
|----------|---------------------------------|
| T01      | ZIR 2.0                         |
| T02      | ZIR 1.0                         |
| T03      | ZIR 0.6                         |
| T04      | 1.5R (implant-hole milling bur) |
| T05      | PMMA 2.0                        |
| T06      | PMMA 1.0                        |
| T07      | PMMA 0.6                        |

### 3.7 Installing and Removing the Fixture

The YRC-5X fixture is an open fixture for material-block installation, as shown below.


Fixture — open position

#### 3.7.1 Installing the Material Block

Seat the  $\varnothing 98$  mm circular block in the fixture as shown.

1. Open the working-cabin door and clean the clamp recess.
2. Use a screwdriver to loosen the three screws on the pressure plate.
3. Put the material block into the fixture and tighten the three fixing screws with a screwdriver, as shown.

Fixture with material block installed

*Note: The screws must be tightened. If the material block moves or vibrates during cutting, it may damage the equipment.*

#### 3.7.2 Removing the Material Block

To remove the block, simply release the three fixing screws and remove it.

## 3.8 Program Interruption and Stop

*Note: When the program is interrupted or stopped, a corresponding alarm window pops up on the operating system page. When air pressure is insufficient, the program automatically pauses and restarts automatically once pressure is sufficient.*

The program must be canceled when:

- Equipment failure occurs.
- The tool is broken.
- A power failure occurs.

If the program is cancelled, it must be restarted on the device.

### 3.8.1 What to Do When the Program Is Interrupted

If an interruption occurs, the operating system page displays a corresponding message. If "insufficient air pressure" is shown, check the following:

- Observe whether the pressure on the pressure-regulating filter element and gauge is within the correct range.
- Whether there is air leakage in the trachea (air hose).
- Whether the compressed-air pump is normal.

### 3.8.2 How to Deal with Equipment Failure

Equipment failures are caused by commands issued by the control unit when an emergency occurs. The operating system page displays error messages and error codes. If the alarm cannot be cleared, record the alarm information promptly and contact after-sales service.

### 3.8.3 How to Deal with Tool Damage

If the tool breaks during processing, the machine cannot recognize this immediately, and the spindle will continue to rotate and cut until the next bur change. Only when the bur length is measured, if the measured length is wrong, will the running program report an error, displayed on the operating system page.

Tool breakage can be caused by:


- The bur is damaged or worn.
- The material selected during typesetting is inconsistent with the actual material being processed.
- The bur is misplaced or manually inserted into the spindle at the wrong time.

In such cases, this bur is not suitable for further processing in this step. Replace it as follows:

1. Open the working chamber.
2. Take out the damaged bur from the working chamber.
3. If the spindle grabbed the bur from the tool magazine, check whether it is placed in the correct position; after confirming, put a new bur in the correct position.

4. If the bur was manually put into the spindle, check whether the damaged bur matches the bur shown on the operating system page; prepare a new, correct bur.
5. Close the working chamber and restart the procedure.

## 4. Introduction to the Operating System Interface

Operating system page — main page

Program selection ("Load") · Download program ("Input") · Frequently used functions ("System Operation") · Coordinate display · Start button.

### 4.1 Status Display


- Current program — the program currently copied into the operating system and being executed.
- Status indicator: green = processing in progress; red = alarm; green & red = processing completed.
- Process time — elapsed running time of the current processing program.
- Current tool number — the tool number currently held by the spindle clamp ("0" if none is clamped).

### 4.2 Button Functions

Page 1 (main page):

1. Breakpoint Process — resume processing after a normal stop during processing.
2. Home — turn on the power and return to the original position.
3. Start-up — after loading is complete, click to have the device automatically run the machining program.

4. **Pause** — pause a running program; click "Start-up" again to resume. Note: "Pause" only temporarily stops the movement of the motion axes; it does not stop spindle rotation.
5. **Reset** — reset all current alarm conditions and restore the device to its original state; also stops the program completely while running. Note: after clicking "Reset", the spindle stops slowly, so wait before opening the door for the next operation. To terminate the equipment immediately in a serious condition, use the Emergency Stop switch.
6. **Load** — select the pending program; all programs in the system are shown, default format is the NC file.
7. **Input** — when copying a program via USB, click "Input" and select the program from the list on the left of the pop-up window to import it.
8. **System Operation** — page-turning function button.

*Operating system page 2 — auxiliary page*

#### **Page 2 (auxiliary page):**

1. **Language selection** (English, Chinese, Russian, Spanish).
2. **A180** — displays the main-axis angle.
3. **Clamp** — click Run to move the fixture spindle to given coordinates to facilitate clamping. Note: entering Clamp, click "P1" to move the spindle to a safe position, rotate axis B to 0°, and set axis A to 0°.
4. **Auto T** — detects the tool length of the bur on the current spindle and records it into the system.
5. **Tool changing** — manually change the tool; automatically completes the tool-change function. Replace burs (T1–T7 selectable).
6. **Tool reset** — preset function; all tools are placed into the magazine and each working axis returns to a safe position.

7. Software update — install the program.
8. Input — import a program.
9. Delete files — when the system has too many stored files, click "Delete files" and select the files to remove in the pop-up window.
10. Network setting — when the device is connected to the network, enter this interface and manually match the correct IP address and gateway to transfer NC files over the network.
11. Blow — click to manually open the spindle blow valve; click again to close.
12. Tool magazine display — shop of milling tools.
13. Spindle On — manually turn on spindle rotation.
14. Blow — manual switch for the spindle air blow.
15. Working lamp indicator light — red indicates a failure alarm state; green indicates normal processing; yellow/green flashing indicates completion of the processing program.
16. Loose/Clamp — spindle three-jaw spring-clip switch.
17. MPG / MPG analog off — functions reserved for manufacturer use only.
18. Maintain — device authority management.
19. << / >> — previous / next page.

## 4.3 Coordinate Display Area

Displays the current mechanical coordinates and the set external deviation coordinates.

- Coordinate — store the coordinates set by the current equipment, or recall previously archived coordinates through this function if the equipment deviates.
- Outer Coordinate Modification ribbon — select the axis to modify, enter the coordinate value, and click "Write in" to overwrite the original coordinate.

## 4.4 Tool Position

Click to open the cutter-base coordinate setting interface.

# 5. Equipment Maintenance and Cleaning

## 5.1 Safety Notes

*Regular maintenance ensures processing quality and prolongs the machine's service life. Make sure each axis is in a safe position, turn off the power, and unplug the power cord before servicing or maintenance.*

**Warning:** *For a short time after power-off, electronic components may still be partially charged — wait a few minutes to ensure full discharge before servicing.*

- Clean the device with a damp cloth and, if necessary, a suitable cleaning agent. Do not use abrasive detergents or agents containing ingredients corrosive to rubber.
- All operations must be carried out by trained, qualified personnel.
- Before working on the pneumatic system, make sure it is turned off and the air pressure in the whole passage has been reduced to normal.

## 5.2 Accessory Replacement

*Warning: Using non-original or unsuitable accessories can put operators at risk and may damage the equipment or render it inoperable. Only use original accessories from the manufacturer, or accessories authorized by the manufacturer. For any questions, contact the manufacturer's customer service.*

## 5.3 Routine Cleaning

- Routine cleaning extends the life of individual parts and avoids failures caused by debris.
- Regular cleaning extends the overall service life of the equipment.
- Clean plastic parts with a suitable liquid cleaner. It is not recommended to use an air gun for interior cleaning, to prevent dust and debris from entering mechanical parts and causing damage.

## 5.4 Material Fixture Cleaning

When placing materials, make sure the clamps and screws are clean to achieve the best grip.

## 5.5 Cleaning the Tool Library and Tools

The tool library and the tools must be free of chips and other debris to ensure accurate length detection and tool replacement. In case of damage, contact the manufacturer immediately.

## 5.6 Cleaning the Spindle Clamp

- Do not spray oil or compressed air directly on the spindle head — this may let these substances penetrate the spindle bearing, affecting rotation or causing bearing damage.
- The spindle chuck and inserted bur must be clean; impurities entering the spindle increase rotation resistance and affect concentricity.
- Normal cleaning cycle of the spindle clamp is weekly. Also clean the spindle clip whenever the pin is damaged by impact.

Spindle chuck hexagon wrench

*Warning: Do not close the spindle chuck when the chuck is not in the spindle — spinning the spindle without a collet can damage it. Always keep a bur on the collet.*

#### Cleaning of the spindle clip:

1. Return the equipment to a safe position.
2. Ensure the spindle is not rotating.
3. Use a special tool to remove the spindle clip.
4. Clean the clip-head mounting hole.
5. Clean the clip with small brushes, ensuring the clip head is free of impurities.
6. Reposition the clean clip back on the spindle, ensuring it is tightened in place.

## 5.7 Checking Hoses, Cables and Connections

1. Power off the equipment and cut off the compressed-air supply.
2. Check the hoses and cables to ensure connections are normal.
3. If there is damage, do not restart the device — replace damaged cables and hoses immediately with spare parts.
4. Make sure connections on the connection board are firm and not loose.
5. Reinsert loose cables and hoses into the appropriate connections. If a connection is damaged, stop using the device immediately and contact after-sales service.

## 5.8 Checking the Pressure-Regulating Filter

*Caution: Contamination of the regulator filter may damage the equipment.*

- Check the pressure-regulating filter for contamination weekly.
- If moisture, oil or other solid particles are found inside, stop using the device.

### 5.8.1 Checking the Regulator Filter for Condensation


Condensation in the pressure-regulating filter usually indicates the compressed air is not dry enough.

1. Check whether water, oil or solid particles have accumulated in the water separator. If so, stop using the machine immediately and proceed to the next step.
2. Check whether the compressed-air supply meets the equipment's operating requirements. Do not use the machine until requirements are met.
3. Pull down the pressure-regulating filter drainer to discharge impurities.

## 5.9 Maintenance of Screw and Guide Rail

Regular maintenance of the screw and guide rail helps maintain machining accuracy and prolongs the equipment's service life.

1. Turn off the power of the device.
2. Remove the fixing screw of the pressure-regulating filter and take off the filter, as shown.
3. Remove the screws of the side cover and rear cover of the device, and remove both covers, as shown.

*Removal of the equipment cover plate*

1. Use a vacuum cleaner to clean the dust inside the equipment; use a clean paper towel to wipe stains off the screw rod and guide rail.
2. Apply grease evenly to the surface of the wiped lead screw and guide rail.
3. Reinstall the removed cover plates and pressure-regulating filter onto the device.

If any problem arises during maintenance, please contact after-sales service promptly.

## 5.10 Maintenance Schedule

| Maintenance items                                        | Cycle / period   |
|----------------------------------------------------------|------------------|
| Clean working cabin, tool holder, tool library, and burs | Every day        |
| Clean the spindle clip                                   | Weekly           |
| Pressure-regulating filter                               | Weekly           |
| Clean screw rods and guide rails                         | Every six months |
| Machine-tool calibration                                 | If necessary     |

## 5.11 Equipment Calibration

The device is calibrated at the factory. As long as the machining result is flawless, no calibration is required. Calibration takes time and, if done incorrectly, can lead — in the worst case — to damage to the device.

If the processing result is not satisfactory, first adjust the working conditions: check the fixation of the workpiece, the state of the bur, and the processing material. Contact the supplier before calibrating the equipment. During calibration, take great care measuring and entering data — when in doubt, stop calibration.

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